

SOP.UCF.027 - IDEAS SOFTWARE TUTORIAL



- 1. Turn on the computer
- Login with your Agendo Credentials:
 Username Your email Password Your Agendo's password
- 3. Open IDEAS Software and click on the Start Analysis option
- On File menu select Open to open your files
 (.rif = Raw image file; .cif = Compensated image file; .daf = Data analysis file)
- 5. Select an analysis template (.ast file)
- 6. Select the used channels and the **Image Gallery Display** properties that will be applied to the images
- 7. If applicable, select the desired Wizard and follow the indicated steps
- To add additional analysis, use the histogram and dot plot icons, or select Guided Analysis, choose Build Blocks, and automatically create plots to the analysis
- 9. When the analysis is complete, save the data file (.daf) and analysis template (.ast)
- 10. When the Wizard is complete, validate the regions by opening the control files using the saved analysis template and compensation matrix. Also verify the data in the statistics report under **Reports** and **Define Statistics Report**

Batch the analysis

- a. Complete you analysis and ensure that everything is correctly analysed
- b. On **Tools** menu, select **Batch Data Files**, then select **Add Batch**, and then **Add Files**. Select all files to be batch processed (you can batch .rif, .cif or .daf)
- c. If data is not already compensated, select the **Compensation Matrix** to be applied
- d. Select the **Template** to apply to the batched data
- e. Choose Submit Batches and wait for the IDEAS to process data

Generate the report

- a. Reporting Images
 - i. Select the Gallery Display icon
 - ii. Select the channel to adjust and optimize the image contrast
 - iii. Create composite overlays
 - iv. Create Gallery view



SOP.UCF.027 - IDEAS SOFTWARE TUTORIAL



- v. Select images for you figure, save them as a population and display in the image gallery
- vi. Right-click and select **Copy / Save Gallery**. Edit display to include measurement tool feature values, object number, channel names and paste into reporting software
- b. Reporting **Plots**
 - i. Right-click on plot and select Graph Properties
 - ii. Select Statistics icon and add / remove statistics
 - iii. Right-click on the plot and select Copy / Save Graph
 - iv. Save to clipboard and paste into reporting software

Histogram overlays

- a. On IDEAS
 - i. Select Tools / merge .cif
 - ii. Crtl-select the data files to overlay
 - iii. Open the merged .cif in the experiment template
 - iv. Create the population to overlay for each sample in the merged data file
 - v. Create a histogram and in the histogram properties Crtl-select each population to overlay
- b. On Microsoft Office Word
 - i. Copy / Paste each histogram into Word
 - ii. Select the top plot and set transparent color and repeat for each histogram making each below visible

Images

- a. To merge **Raw Image Files**, on Tools menu, select **Merge .rif Files** and **add files** to merge
- b. To save an **analysis template** file, from the File drop down option select Save as template file (.ast)



SOP.UCF.027 - IDEAS SOFTWARE TUTORIAL



Compensation

- a. To create a compensation matrix
 - i. Select New matrix when opening a .rif select Compensation, then Create New Matrix
 - Select the control files for compensation, which are the ones that are record as "<u>-noBF</u>". Click **next** to load the files
 - iii. Validate the compensation matrix, checking for high values, that will be highlighted in red and optimize the matrix, correcting first the biggest ones:
 - I. Double-click the red value in the matrix and the Matrix coefficient intensity plot is displayed
 - II. Use a region tool to select a new positive population that excludes any outlier and assign the new population to its channel
 - III. A new matrix is calculated and should be re-evaluated
- b. Click **Finish** and save the compensation matrix